

# Use Cases

## 1. Use Case: CES Device Session

**Primary Actor:** CES Device User

**Stakeholders and Interests:**

- CES Device User - wants to use Cranio-Electro stimulation to treat their condition(s).
- Manufacturing Company - wants to make sure that the electrotherapy is working as intended.

**Precondition:**

There is enough battery to activate the CES device. The device must also be powered on. The connection level should be sufficient (Okay, Good).

**Success Guarantee:**

The device emits an electrotherapeutic treatment for the desired length of time, or until the battery is no longer charged.

**Main Success Scenario:**

1. Turn on device (include Power On)
2. CES Device battery level is shown with a battery percentage of 100%.
3. User presses the power button again to select the time period (twenty-minutes, forty-five-minutes).
4. User choose whether to record session (include Save Therapy)
5. CES does a routine check (include Check Connection)
6. User chooses a session (ALPHA, THETA, DELTA, MET).
7. Session is confirmed with the start session button.
8. User adjusts intensity until the stimulation can barely be felt (include Change Intensity).

**Extensions**

3a. Arbitrary session length.

- The CES device user can select the length in minutes button, and assign an arbitrary amount of time for a custom session.

## 2. Use Case: Power On

**Primary Actor:** CES Device User

**Stakeholders and Interests:**

- CES Device User - wants to use Cranio-Electro stimulation to treat their condition(s).
- Manufacturing Company - wants to make sure that the electrotherapy is working as intended.

**Precondition:**

CES Device user has a device with sufficient battery power.

**Success Guarantee:**

The device will startup into the main GUI.

**Main Success Scenario:**

1. User pushes the power button to turn on the CES device.
2. The device turns on.

**Extensions:**

- 1a. Battery level is low.
  - Battery low warnings will be shown visually as a function of the length and intensity.

## 3. Use Case: Power Off

**Primary Actor:** CES Device User

**Stakeholders and Interests:**

- CES Device User - wants to use Cranio-Electro stimulation to treat their condition(s).
- Manufacturing Company - wants to make sure that the electrotherapy is working as intended.

**Precondition:**

CES Device user has a device that is powered on.

**Success Guarantee:**

The device will shut down.

**Main Success Scenario:**

1. User pushes the power button to turn off the CES device.
2. The device turns off.

## 4. Use Case: Save Therapy

**Primary Actor:** CES Device User

**Stakeholders and Interests:**

- CES Device User - wants to use Cranio-Electro stimulation to treat their condition(s).
- Manufacturing Company - wants to make sure that the electrotherapy is working as intended.

**Precondition:**

CES Device user has a device that is powered on. Adequate connection level (Okay, Good).

**Success Guarantee:**

The record is appended to the SavedTherapies.txt logging file.

**Main Success Scenario:**

1. User clicks on the “Record Session” radio button.
2. After the treatment, the user can add the therapy to persistent memory, containing the session type, duration, and intensity level.

**Extensions:**

- 2a. The CES device session finishes.
  - After a completed session the therapy information is saved in the logging file.
- 2b. The CES device runs out of battery.
  - After an abrupt power off the device saves the record, before shutting off.

## 5. Use Case: Change Intensity

**Primary Actor:** CES Device User

**Stakeholders and Interests:**

- CES Device User - wants to use Cranio-Electro stimulation to treat their condition(s).
- Manufacturing Company - wants to make sure that the electrotherapy is working as intended.

**Precondition:**

CES Device user has a device that is powered on. There is sufficient power to change to the specified intensity. Adequate connection level (Okay, Good).

**Success Guarantee:**

The intensity meter changes to the specified level.

**Main Success Scenario:**

1. User presses the decrement/increment button to adjust the intensity of the CES device.
2. The UI is updated to represent the intensity change.

## 6. Use Case: Check Connection

**Primary Actor:** CES Device User

**Stakeholders and Interests:**

- CES Device User - wants to use Cranio-Electro stimulation to treat their condition(s).
- Manufacturing Company - wants to make sure that the electrotherapy is working as intended.

**Precondition:**

CES Device user has a device that is powered on. There is sufficient power to change to the specified intensity.

**Success Guarantee:**

The intensity meter changes to the specified level.

**Main Success Scenario:**

1. CES Device does a connection test.

**Extension:**

1a. *“Left Ear Disconnected”*:

- The left ear is disconnected, causing the session to stop momentarily until reconnected.

1b. *“Right Ear Disconnected”*:

- The right ear is disconnected, causing the session to stop momentarily until reconnected.

1c. *“Left and Right Ear Disconnected”*:

- The left and right ear is disconnected, causing the session to stop momentarily until reconnected.

1d. *“No Connection”*:

- If there is no connection to the device for 20 seconds, the device stops stalling and powers off.

1e. *“Okay Connection”*:

- If the connection is *“Okay”* then the session proceeds, though drains more battery than if the connection was *“Excellent”*.

1e. *“Excellent Connection”*:

- If the connection is *“Excellent”* then the device proceeds with minimal entropic loss.

## Use Case Diagram: CES Device

