

Use Case: Power On

Primary Actor: CES Device User

Stakeholders and Interests:

- CES Device User - wants to use Cranio-Electro stimulation to treat their condition(s).
- Manufacturing Company - wants to make sure that the electrotherapy is working as intended.

Precondition:

CES Device user has a device with sufficient battery power.

Success Guarantee:

The device will startup into the main GUI.

Main Success Scenario:

1. User pushes the power button to turn on the CES device.
2. The device turns on.

Extensions:

- 1a. Battery level is low.
 - Battery low warnings will be shown visually as a function of the length and intensity.