

Use Case: Power Off

Primary Actor: CES Device User

Stakeholders and Interests:

- CES Device User - wants to use Cranio-Electro stimulation to treat their condition(s).
- Manufacturing Company - wants to make sure that the electrotherapy is working as intended.

Precondition:

CES Device user has a device that is powered on.

Success Guarantee:

The device will shut down.

Main Success Scenario:

1. User pushes the power button to turn off the CES device.
2. The device turns off.

Extensions:

2a. If the session is being recorded and the device has to power off, the CES device will save the session before exiting.